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
1. Maldonado-Codina C, Morgan PB. In vitro water wettability of silicone hydrogel contact lenses determined using the sessile drop and captive bubble techniques. *J Biomed Mater Res A*. 2007; 83 (2): 496–502.
2. Jones L, Subbaraman L, Rogers R, Dumbleton K. Surface treatment, wetting and modulus of silicone hydrogels. *Optician*. 2006; 232 (6067): 28–34.
3. Weikart CM, Matsuzawa Y, Winterton L, Yasuda HK. Evaluation of plasma polymer-coated contact lenses by electrochemical impedance spectroscopy. *J Biomed Mater Res* 2001; 54 (4): 597–607.
4. Gromacki SJ. Caring for silicone hydrogel contact lenses – Part 2. *Contact Lens Spectrum*; June 2005.
5. Lopez-Aleman A, Compan V, Refojo MF. Porous structure of Purevision versus Focus Night & Day and conventional hydrogel contact lenses. *J Biomed Mater Res* 2002; 63 (3): 319–325.
6. Maiden AC, Vanderlaan DG, Turner DC, Love RN, Ford JD, Molock FF, Steffen RB, Hill GA, Alli A, McCabe KP, inventors; Johnson & Johnson Vision Care, Inc., assignee. Biomedical devices containing internal wetting agents. United States patent US 6,367,929. April 9, 2002.
7. Wygladacz K, Hook D, Steffen R, Reindel W. Breaking the cycle of discomfort. *Contact Lens Spectrum*; special edition 2014; 29 (13): 23–28.
8. Jones L, May C, Nazar L, Simpson T. In vitro evaluation of the dehydration characteristics of silicone hydrogel and conventional hydrogel contact lens materials. *Cont Lens Anterior Eye*. 2002; 25 (3): 147–156.
9. Bhamla MS, Nash WL, Elliott S, Fuller GG. Influence of lipid coatings on surface wettability characteristics of silicone hydrogels. *Langmuir*. 2015; 31 (13): 3820–3828.
10. Lin MC, Svitova TF. Contact lenses wettability in vitro: effect of surface-active ingredients. *Optom Vis Sci*. 2010; 87 (6): 440–447.
11. Tonge S, Jones L, Goodall S, Tighe B. The ex vivo wettability of soft contact lenses. *Curr Eye Res*. 2001; 23 (1): 51–59.
12. Ketelson HA, Meadows DM, McQueen N, Stone RP. Enhancing wettability with multi-purpose solutions. *Rev Cornea Cont Lens*. 2005; 142 (1 Suppl): 44–47.
13. Aragona P. Hyaluronan in the treatment of ocular surface disorders (Chapter 24). In: HG Garg and CA Hales, eds. *Chemistry and Biology of Hyaluronan*. Oxford: Elsevier Ltd; 2004 (pages 529–551).
14. Necas J, Bartosikova L, Brauner P, Kolar J. Hyaluronic acid (hyaluronan): a review. *Vet Med* 2008; 53: 397–411.
15. Brown TJ, Alcorn D, Fraser JR. Absorption of hyaluronan applied to the surface of intact skin. *J Invest Dermatol*. 1999; 113 (5): 740–746.
16. Johnston JB. A simple, nondestructive assay for bound hyaluronan. *J Biomed Mater Res*. 2000; 53 (2): 188–191.
17. Laurent TC, Laurent UB, Fraser JR. Functions of hyaluronan. *Ann Rheum Dis*. 1995; 54 (5): 429–432.
18. FDA 510 (k) Summary K121201. Soft (Hydrophilic) Contact Lens, Oct 26, 2012.
19. FDA 510 (k) Summary K083757. Bausch & Lomb BP-ZO2 Multi-Purpose Solution. November 18, 2009.
20. FDA 510k Summary K032030. blink™ CL Lubricant Eye Drops, Sep. 25, 2003.
21. Lerner LE, Schwartz DM, Hwang DG, Howes EL, Stern R. Hyaluronan and CD44 in the human cornea and limbal conjunctiva. *Exp Eye Res*. 1998; 67 (4): 481–484.
22. Wysesbeek YS, Loya N, Ben Sira I, Ophir I, Ben Shaul Y. The effect of sodium hyaluronate on the corneal epithelium. An ultrastructural study. *Invest Ophthalmol Vis Sci*. 1988; 29 (2): 194–199.

23. Johnson ME, Murphy PJ, Boulton M. Effectiveness of sodium hyaluronate eyedrops in the treatment of dry eye. *Graef Arch Clin Exp Ophthalmol*. 2006; 244 (1): 109–112.
24. Mengher LS, Pandher KS, Bron AJ, Davey CC. Effect of sodium hyaluronate (0.1%) on break-up time (NIBUT) in patients with dry eyes. *Br J Ophthalmol*. 1986; 70 (6): 442–447.
25. McCann LC, Tomlinson A, Pearce EI, Papa V. Effectiveness of artificial tears in the management of evaporative dry eye. *Cornea*. 2012; 31 (1): 1–5.
26. Condon PI, McEwen CG, Wright M, Mackintosh G, Prescott RJ, McDonald C. Double blind, randomised, placebo controlled, crossover, multicenter study to determine the efficacy of a 0.1% (w/v) sodium hyaluronate solution (Fermavisc) in the treatment of dry eye syndrome. *Br J Ophthalmol*. 1999; 83 (10): 1121–1124.
27. Gonzalez-Mejome JM, Lopez-Aleman A, Almeida JB, Parafita MA, Refojo MF. Microscopic observation of unworn siloxane-hydrogel soft contact lenses by atomic force microscopy. *J Biomed Mater Res B Appl Biomater*. 2006; 76 (2): 412–418.
28. Teichroeb JH, Forrest JA, Ngai V, Martin JW, Jones L, Medley J. Imaging protein deposits on contact lens materials. *Optom Vis Sci*. 2008; 85 (12): 1151–1164.
29. Scheuer CA, Fridman KM, Barniak VL, Burke SE, Venkatesh S. Retention of conditioning agent hyaluronan on hydrogel contact lenses. *Cont Lens Anterior Eye*. 2010; 33 (Suppl 1): S2–S6.
30. Liu XM, Harmon PS, Maziarz EP, Rah MJ, Merchea MM. Comparative studies of hyaluronan in marketed ophthalmic products. *Optom Vis Sci*. 2014; 91 (1): 32–38.
31. Scheuer C, Rah MJ, Reindel WT. Increased concentration of hyaluronan in tears after soaking contact lens in Biotrue multipurpose solution. *Clin Ophthalmol*. 2016, in press.
32. Tas J. The Alcian blue and combined Alcian blue-Safranin O staining of glycosaminoglycans studied in a model system and in mast cells. *Histochem J*. 1977; 9 (2): 205–230.
33. Luo Y, Prestwich GD. Hyaluronic acid-N-hydroxysuccinimide: a useful intermediate for bioconjugation. *Bioconjug Chem*. 2001; 12 (6): 1085–1088.
34. Homer KA, Denbow L, Whiley RA, Beighton D. Chondroitin sulfate depolymerase and hyaluronidase activities of viridans streptococci determined by a sensitive spectrophotometric assay. *J Clin Microbiol*. 1993; 31 (6): 1648–1651.
35. Fagnola M, Pagani MP, Maffioletti S, Tavazzi S, Papagni A. Hyaluronic acid in hydrophilic contact lenses: spectroscopic investigation of the content and release in solution. *Cont Lens Anterior Eye*. 2009; 32 (3): 108–112.
36. Tavazzi S, Tonveronachi M, Fagnola M, et al. Wear effects on microscopic morphology and hyaluronan uptake in siloxane-hydrogel contact lenses. *J Biomed Mater Res B Appl Biomater*. 2015; 103 (5): 1092–1098.
37. Wygladacz K, Bakker E. Fluorescent microsphere fiber optic micro-sensor array for direct iodide detection at low picomolar concentrations. *Analyst*. 2007; 132 (3): 268–272.
38. Khan AY, Saha B, Suresh Kumar G. Phenazinium dyes safranin O and phenosafranin induce self-structure in single stranded polyadenylic acid: structural and thermodynamic studies. *J Photochem Photobiol B*. 2014; 132: 17–26.
39. Weeks A, Boone A, Luensmann D, Jones L, Sheardown H. The effects of hyaluronic acid incorporated as a wetting agent on lysozyme denaturation in model contact lens materials. *J Biomater Appl*. 2013; 28 (3): 323–333.
40. Singh A, Li P, Beachley V, McDonnell P, Elisseeff JH. A hyaluronic acid-binding contact lens with enhanced water retention. *Cont Lens Anterior Eye*. 2015; 38 (2): 79–84.

Refreshing the contact lens surface

Biotrue multipurpose solution (MPS) includes hyaluronic acid (HA) as a natural wetting agent. In an in vitro study, the morphology and distribution of HA sorbed from both Biotrue and pure HA solution of equal concentration (0.1%) on SiHy contact lens surfaces was examined. Atomic force microscopy imaging was used to illustrate the topography of SiHy contact lens materials. The distribution, as well as fine details of the HA network, were resolved by staining HA with Gram's safranin, then imaging with confocal laser-scanning microscopy and differential interference-contrast microscopy. In this approach, SiHy materials from different generations was chosen for the study: Balafilcon A with plasma treatment of surface, Senofilcon A and Samfilcon A as SiHy materials fabricated with an internal wetting agent (polyvinylpyrrolidone – PVP). Studies demonstrated that HA sorbed from Biotrue MPS on different types of silicone hydrogel (SiHy) contact lens materials.

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